

Protein extraction from *Nicotiana benthamiana* leaf Institute of Genetics, Heinrich-Heine University Düsseldorf, Germany

CONTEXT

The laboratory is focused on developmental biology of plants using the model system *Arabidopsis thaliana*.

The aim of this study was to investigate the intracellular localization of protein receptors (CLV1, CLV2, and CRN) in plant cells and their tendencies for protein-protein interactions. To analyze receptor localization and interaction, a transient expression system in *Nicotiana benthamiana* leaf epidermis cells was used [1].

MATERIAL

- Precellys®24 homogenizer.
- Precellys® kit: 03961-1-003 (ceramic beads 1.4mm)
- Sample : ~0.1 g of *N. benthamiana* leaf tissue.
- Extraction buffer: 750 µL (0.1 M Tris-HCl, pH 8.3, 5 mM dithiothreitol, 5 mM EDTA and protease inhibitor).

PROTOCOL

- Plant tissue homogenized using the Precellys®24: 5500 rpm, 1x20 sec.
- 1h incubation of plant extract at 4°C followed by 10 min denaturation at 95°C.
- Protein separation by SDS-PAGE & Western-blot analysis.

[1] A. Bleckmann et al., Plant Physiology, January 2010, Vol. 152, pp. 166–176, www.plantphysiol.org



RESULTS

Our transient expression studies of translational fusions with GFP or mCherry now showed that all three receptor proteins can localize to the plasma membrane and have the capacity to undergo multiple interactions.




Figure 1: Western-blot analysis of protein extracts from *N. benthamiana* leaf cells transiently expressing CLV1-GFP, CLV2-GFP, or CRN-GFP. An anti-GFP antibody was used for detection; sizes of protein markers are given in kD. The Ponceau S-stained protein bands of Rubisco are shown as a loading control. wt, wild type.

The Precellys®24 is a fast method for protein isolation. Protein extracts were used to demonstrate fusion protein expression and stability.

CONCLUSION

Precellys®24 is well adapted for homogenization of plant tissues to extract and study proteins expression. Sample preparation is not only easy but cross contamination free.

The Precellys®24 can also be used for efficient isolation of DNA and RNA from plant tissue.

For more details, please contact
precellys@bertin.fr

